

CLVS2 Antibody, Biotin conjugated



PACO36325

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Required for normal morphology of late endosomes and/or lysosomes in neurons By similarity. Binds phosphatidylinositol 3,5-bisphosphate (PtdIns(3,5)P2).

Gene ID:

CLVS2

Uniprot

Q5SYC1

Synonyms:

Clavesin-2 (Retinaldehyde-binding protein 1-like 2) (clathrin vesicle-associated Sec14 protein 2), CLVS2, C6orf212 C6orf213 RLBP1L2

Immunogen:

Recombinant Human Clavesin-2 protein (1-327AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A