

GLT8D1 Antibody



PACO36342

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:

Golgi apparatus, membrane, transferase activity, transferring glycosyl groups, polysaccharide biosynthetic process.

Gene ID:

GLT8D1

Uniprot

Q68CQ7

Synonyms:

Glycosyltransferase 8 domain-containing protein 1 (EC 2.4.1. -), GLT8D1, GALA4A


Immunogen:

Recombinant Human Glycosyltransferase 8 domain-containing protein 1 protein (33-332AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: GLT8D1 antibody at 10 μ g/ml + Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 42, 29 kDa. Observed band size: 42 kDa.



Immunohistochemistry of paraffin-embedded human liver tissue using PACO36342 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36342 at dilution of 1:100.