

NAIF1 Antibody, HRP conjugated



PACO36347

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Induces apoptosis.

Gene ID:

NAIF1

Uniprot

Q69YI7

Synonyms:

Nuclear apoptosis-inducing factor 1, NAIF1, C9orf90

Immunogen:

Recombinant Human Nuclear apoptosis-inducing factor 1 protein (1-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A