

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Polyglutamylase which preferentially modifies alpha-tubulin. Involved in the side-chain initiation step of the polyglutamylation reaction rather than in the elongation step. Required for CCSAP localization to both spindle and cilia microtubules. Increases the effects of NCOA2 in glucocorticoid receptor-mediated repression and induction and in androgen receptor-mediated induction.

**Gene ID:**

TTLL5

**Uniprot**

Q6EMB2

**Synonyms:**

Tubulin polyglutamylase TTLL5 (EC 6. -. -. -) (SRC1 and TIF2-associated modulatory protein) (Tubulin--tyrosine ligase-like protein 5), TTLL5, KIAA0998 STAMP

**Immunogen:**

Recombinant Human Tubulin polyglutamylase TTLL5 protein (1-269AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A