

C12orf60 Antibody



PACO36414

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:**Gene ID:**

C12orf60

Uniprot

Q5U649

Synonyms:

Uncharacterized protein C12orf60, C12orf60


Immunogen:

Recombinant Human Uncharacterized protein C12orf60 protein (1-245AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36414 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO36414 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO36414 at dilution of 1:100.