

C10orf91 Antibody



PACO36450

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C10orf91

Uniprot

Q5T1B1

Synonyms:

Uncharacterized protein C10orf91, C10orf91


Immunogen:

Recombinant Human Uncharacterized protein C10orf91 protein (1-145AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO36450 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO36450 at dilution of 1:100.