

EGFLAM Antibody, Biotin conjugated



PACO36469

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in both the retinal photoreceptor ribbon synapse formation and physiological functions of visual perception. Necessary for proper bipolar dendritic tip apposition to the photoreceptor ribbon synapse. Promotes matrix assembly and cell adhesiveness.

Gene ID:

EGFLAM

Uniprot

Q63HQ2

Synonyms:

Pikachurin (Agrin-like protein) (EGF-like, fibronectin type-III and laminin G-like domain-containing protein), EGFLAM, AGRINL AGRNL PIKA

Immunogen:

Recombinant Human Pikachurin protein (25-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A