
Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Putative nucleoside deaminase. May catalyze the hydrolytic deamination of adenosine or some similar substrate and play a role in purine metabolism.

Gene ID:

ADAL

Uniprot

Q6DHF7

Synonyms:

Adenosine deaminase-like protein (EC 3.5.4. -), ADAL


Immunogen:

Recombinant Human Adenosine deaminase-like protein (1-267AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: ADAL antibody at 6 μ g/ml + Mouse kidney tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 41, 31, 38 kDa. Observed band size: 50 kDa.