

# ADAL Antibody, HRP conjugated



PACO36483

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Putative nucleoside deaminase. May catalyze the hydrolytic deamination of adenosine or some similar substrate and play a role in purine metabolism.

**Gene ID:**

ADAL

**Uniprot**

Q6DHV7

**Synonyms:**

Adenosine deaminase-like protein (EC 3.5.4. -), ADAL

**Immunogen:**

Recombinant Human Adenosine deaminase-like protein (1-267AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A