

ACBD5 Antibody, FITC conjugated



PACO36492

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Acyl-CoA binding protein which acts as the peroxisome receptor for pexophagy but is dispensable for aggrephagy and nonselective autophagy. Binds medium- and long-chain acyl-CoA esters.

Gene ID:

ACBD5

Uniprot

Q5T8D3

Synonyms:

Acyl-CoA-binding domain-containing protein 5, ACBD5, KIAA1996

Immunogen:

Recombinant Human Acyl-CoA-binding domain-containing protein 5 protein (1-461AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A