

# ACBD5 Antibody, Biotin conjugated



PACO36493

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Acyl-CoA binding protein which acts as the peroxisome receptor for pexophagy but is dispensable for aggrephagy and nonselective autophagy. Binds medium- and long-chain acyl-CoA esters.

**Gene ID:**

ACBD5

**Uniprot**

Q5T8D3

**Synonyms:**

Acyl-CoA-binding domain-containing protein 5, ACBD5, KIAA1996

**Immunogen:**

Recombinant Human Acyl-CoA-binding domain-containing protein 5 protein (1-461AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A