

FAM160B1 Antibody



PACO36502

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

FAM160B1

Uniprot

Q5W0V3

Synonyms:

Protein FAM160B1, FAM160B1, KIAA1600


Immunogen:

Recombinant Human Protein FAM160B1 protein (1-738AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO36502 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon tissue using PACO36502 at dilution of 1:100.