

ORF2 Antibody



PACO36510

Product Information

Size:

50ug

Reactivity:

Hepatitis E virus genotype 1

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Isoform Secreted protein ORF2: Plays a role in the inhibition of host antibody-mediated neutralization without blocking viral cell entry.

Gene ID:

ORF2

Uniprot

Q68985

Synonyms:

Secreted protein ORF2 (Protein ORF2) (pORF2) ORF2

Immunogen:

Recombinant Hepatitis E virus genotype 1 Secreted protein ORF2 protein (20-660AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A