

# DMRTA1 Antibody, FITC conjugated



PACO36560

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

nucleus, core promoter proximal region sequence-specific DNA binding, protein homodimerization activity, transcription factor activity, sequence-specific DNA binding, sex differentiation.

**Gene ID:**

DMRTA1

**Uniprot**

Q5VZB9

**Synonyms:**

Doublesex- and mab-3-related transcription factor A1, DMRTA1, DMO

**Immunogen:**

Recombinant Human Doublesex- and mab-3-related transcription factor A1 protein (201-504AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A