

OLFML2A Antibody



PACO36582

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

OLFML2A

Uniprot

Q68BL7

Synonyms:

Olfactomedin-like protein 2A (Photomedin-1), OLFML2A


Immunogen:

Recombinant Human Olfactomedin-like protein 2A protein (1-300AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO36582 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36582 at dilution of 1:100.