

AVPI1 Antibody, FITC conjugated



PACO36592

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in MAP kinase activation, epithelial sodium channel (ENaC) down-regulation and cell cycling.

Gene ID:

AVPI1

Uniprot

Q5T686

Synonyms:

Arginine vasopressin-induced protein 1 (AVP-induced protein 1), AVPI1

Immunogen:

Recombinant Human Arginine vasopressin-induced protein 1 protein (1-147AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A