

# Trim21 Antibody



PACO36614

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

Trim21

**Uniprot**

Q62191

**Synonyms:**

E3 ubiquitin-protein ligase TRIM21 (EC 2.3.2.27) (52 kDa Ro protein) (52 kDa ribonucleoprotein autoantigen Ro/SS-A) (RING-type E3 ubiquitin transferase TRIM21) (Ro(SS-A)) (Sjogren syndrome type A antigen) (SS-A) (Tripartite motif-containing protein 21), Trim21, Ro52 Ssa1

**Immunogen:**

Recombinant Mouse E3 ubiquitin-protein ligase TRIM21 protein (1-470AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A