

Idh1 Antibody, Biotin conjugated



PACO36633

Product Information

Size:

50ug

Reactivity:

Lactobacillus plantarum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

Idh1

Uniprot

P56512

Synonyms:

L-lactate dehydrogenase 1 (L-LDH 1) (EC 1.1.1.27), Idh1, I-IdhL Idh IdhL IdhL1

Immunogen:

Recombinant Lactobacillus plantarum L-lactate dehydrogenase 1 protein (1-320AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A