

FNDC7 Antibody



PACO36654

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FNDC7

Uniprot

Q5VTL7

Synonyms:

Fibronectin type III domain-containing protein 7, FNDC7

Immunogen:

Recombinant Human Fibronectin type III domain-containing protein 7 protein (234-533AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A