

EFCAB12 Antibody



PACO36694

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

EFCAB12

Uniprot

Q6NXP0

Synonyms:

EF-hand calcium-binding domain-containing protein 12, EFCAB12, C3orf25


Immunogen:

Recombinant Human EF-hand calcium-binding domain-containing protein 12 protein (1-200AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36694 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO36694 at dilution of 1:100.