

OPN5 Antibody, FITC conjugated



PACO36700

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

OPN5

Uniprot

Q6U736

Synonyms:

Opsin-5 (G-protein coupled receptor 136) (G-protein coupled receptor PGR12) (Neuropsin) (Transmembrane protein 13), OPN5, GPR136 PGR12 TMEM13

Immunogen:

Recombinant Human Opsin-5 protein (310-353AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A