

SFXN4 Antibody



PACO36758

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200, IF:1:50-1:200

Protein Background:

Potential iron transporter.

Gene ID:

SFXN4

Uniprot

Q6P4A7

Synonyms:

Sideroflexin-4 (Breast cancer resistance marker 1), SFXN4, BCMR1


Immunogen:

Recombinant Human Sideroflexin-4 protein (2-110AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Western blot. All lanes: SFXN4 antibody at 0.3 μ g/ml. Lane 1: K562 whole cell lysate. Lane 2: U251 whole cell lysate. Lane 3: U87 whole cell lysate. Lane 4: A549 whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 38, 35, 25 kDa. Observed band size: 38 kDa.



Immunohistochemistry of paraffin-embedded human thymus tissue using PACO36758 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO36758 at dilution of 1:100.