

PF14_0076 Antibody, FITC conjugated



PACO36800

Product Information

Size:

50ug

Reactivity:

Plasmodium falciparum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Participates in the digestion of the host hemoglobin. Initial cleavage at the hinge region of hemoglobin, than cleaves at other sites, leading to denaturation of the molecule and to further degradation.

Gene ID:

PF14_0076

Uniprot

Q7KQM4

Synonyms:

Plasmepsin-1 (EC 3.4.23.38) (Aspartic hemoglobinase I) (PfAPG) PF14_0076

Immunogen:

Recombinant Plasmodium falciparum Plasmepsin-1 protein (125-452AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A