

TMPRSS7 Antibody



PACO36810

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Serine protease which preferentially hydrolyzes peptides with Arg at the P1 position.

Gene ID:

TMPRSS7

Uniprot

Q7RTY8

Synonyms:

Transmembrane protease serine 7 (EC 3.4.21. -) (Matriptase-3), TMPRSS7


Immunogen:

Recombinant Human Transmembrane protease serine 7 protein (101-400AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO36810 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO36810 at dilution of 1:100.