

TMPRSS7 Antibody, HRP conjugated



PACO36811

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Serine protease which preferentially hydrolyzes peptides with Arg at the P1 position.

Gene ID:

TMPRSS7

Uniprot

Q7RTY8

Synonyms:

Transmembrane protease serine 7 (EC 3.4.21. -) (Matriptase-3), TMPRSS7

Immunogen:

Recombinant Human Transmembrane protease serine 7 protein (101-400AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A