

LRRN1 Antibody



PACO36854

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

LRRN1

Uniprot

Q6UXK5

Synonyms:

Leucine-rich repeat neuronal protein 1 (Neuronal leucine-rich repeat protein 1) (NLRR-1), LRRN1, KIAA1497


Immunogen:

Recombinant Human Leucine-rich repeat neuronal protein 1 protein (26-631AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO36854 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO36854 at dilution of 1:100.