

Defb14 Antibody



PACO36866

Product Information

Size:

50ug

Reactivity:

Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Has antibacterial activity.

Gene ID:

Defb14

Uniprot

Q7TNV9

Synonyms:

Beta-defensin 14 (BD-14) (mBD-14) (Defensin, beta 14), Defb14


Immunogen:

Recombinant Mouse Beta-defensin 14 protein (23-67AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: Defb14 antibody at 12 μ g/ml + Rat liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 8 kDa. Observed band size: 8 kDa.