

# SNRNP48 Antibody



PACO36874

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Likely involved in U12-type 5' splice site recognition.

**Gene ID:**

SNRNP48

**Uniprot**

Q6IEG0

**Synonyms:**

U11/U12 small nuclear ribonucleoprotein 48 kDa protein (U11/U12 snRNP 48 kDa protein) (U11/U12-48K), SNRNP48, C6orf151

**Immunogen:**

Recombinant Human U11/U12 small nuclear ribonucleoprotein 48 kDa protein (1-300AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A