

ARHGAP11A Antibody



PACO36890

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, regulation of small GTPase mediated signal transduction, small GTPase mediated signal transduction.

Gene ID:

ARHGAP11A

Uniprot

Q6P4F7

Synonyms:

Rho GTPase-activating protein 11A (Rho-type GTPase-activating protein 11A), ARHGAP11A, KIAA0013

Immunogen:

Recombinant Human Rho GTPase-activating protein 11A protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A