

# GFRAL Antibody, HRP conjugated



PACO36907

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Brainstem-restricted receptor for GDF15 which regulates food intake, energy expenditure and body weight in response to metabolic and toxin-induced stresses. Upon interaction with its ligand, GDF15, interacts with RET and induces cellular signaling through activation of MAPK- and AKT- signaling pathways.

**Gene ID:**

GFRAL

**Uniprot**

Q6UXV0

**Synonyms:**

GDNF family receptor alpha-like, GFRAL, C6orf144

**Immunogen:**

Recombinant Human GDNF family receptor alpha-like protein (19-351AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A