

GFRAL Antibody, Biotin conjugated



PACO36909

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Brainstem-restricted receptor for GDF15 which regulates food intake, energy expenditure and body weight in response to metabolic and toxin-induced stresses. Upon interaction with its ligand, GDF15, interacts with RET and induces cellular signaling through activation of MAPK- and AKT- signaling pathways.

Gene ID:

GFRAL

Uniprot

Q6UXV0

Synonyms:

GDNF family receptor alpha-like, GFRAL, C6orf144

Immunogen:

Recombinant Human GDNF family receptor alpha-like protein (19-351AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A