

DHX30 Antibody, HRP conjugated



PACO36923

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an important role in the assembly of the mitochondrial large ribosomal subunit. Required for optimal function of the zinc-finger antiviral protein ZC3HAV1. Associates with mitochondrial DNA.

Gene ID:

DHX30

Uniprot

Q7L2E3

Synonyms:

Putative ATP-dependent RNA helicase DHX30 (EC 3.6.4.13) (DEAH box protein 30), DHX30, DDX30 KIAA0890

Immunogen:

Recombinant Human Putative ATP-dependent RNA helicase DHX30 protein (694-1000AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A