

# NANOGP8 Antibody, HRP conjugated



PACO36959

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May act as a transcription regulator. When overexpressed, promotes entry of cells into S phase and cell proliferation.

**Gene ID:**

NANOGP8

**Uniprot**

Q6NSW7

**Synonyms:**

Homeobox protein NANOGP8, NANOGP8

**Immunogen:**

Recombinant Human Homeobox protein NANOGP8 protein (1-305AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A