

KIRREL2 Antibody



PACO36982

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

May regulate basal insulin secretion.

Gene ID:

KIRREL2

Uniprot

Q6UWL6

Synonyms:

Kin of IRRE-like protein 2 (Kin of irregular chiasm-like protein 2) (Nephrin-like protein 3), KIRREL2, NEPH3


Immunogen:

Recombinant Human Kin of IRRE-like protein 2 protein (21-219AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: KIRREL2 antibody at 8 μ g/ml. Lane 1: Mouse gonadal tissue. Lane 2: Mouse lung tissue. Lane 3: Mouse heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 76, 68, 62, 24 kDa. Observed band size: 76 kDa.