

# BTN2A1 Antibody, Biotin conjugated



PACO36997

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

extracellular exosome, integral component of plasma membrane, lipid metabolic process.

**Gene ID:**

BTN2A1

**Uniprot**

Q7KYR7

**Synonyms:**

Butyrophilin subfamily 2 member A1, BTN2A1, BT2.1 BTF1

**Immunogen:**

Recombinant Human Butyrophilin subfamily 2 member A1 protein (29-248AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A