

FAM171B Antibody



PACO37010

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Gene ID:

FAM171B

Uniprot

Q6P995

Synonyms:

Protein FAM171B, FAM171B, KIAA1946


Immunogen:

Recombinant Human Protein FAM171B protein (59-353AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4


Product Images



Immunohistochemistry analysis of human melanoma using PACO37010 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO37010 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry analysis of human bladder cancer using PACO37010 at dilution of 1:100.