

inIA Antibody



PACO37038

Product Information

Size:

50ug

Reactivity:

Listeria monocytogenes serotype 4b

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Mediates the entry of Listeria monocytogenes into cells.

Gene ID:

inIA

Uniprot

Q723K6

Synonyms:

Internalin-A, inIA


Immunogen:

Recombinant Listeria monocytogenes serotype 4b Internalin-A protein (562-770AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: inIA antibody at 10 μ g/ml + 293T whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 87 kDa. Observed band size: 87 kDa.