

LYPD2 Antibody



PACO37094

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

LYPD2

Uniprot

Q6UXB3

Synonyms:

Ly6/PLAUR domain-containing protein 2, LYPD2, LYPDC2


Immunogen:

Recombinant Human Ly6/PLAUR domain-containing protein 2 protein (23-103AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37094 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO37094 at dilution of 1:100.