

C8orf74 Antibody



PACO37142

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

C8orf74

Uniprot

Q6P047

Synonyms:

Uncharacterized protein C8orf74, C8orf74


Immunogen:

Recombinant Human Uncharacterized protein C8orf74 protein (1-294AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human heart tissue using PACO37142 at dilution of 1:100.