

C16orf87 Antibody, Biotin conjugated



PACO37161

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

C16orf87

Uniprot

Q6PH81

Synonyms:

UPF0547 protein C16orf87, C16orf87

Immunogen:

Recombinant Human UPF0547 protein C16orf87 protein (1-154AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A