

TRAT1 Antibody



PACO37162

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Stabilizes the TCR (T-cell antigen receptor)/CD3 complex at the surface of T-cells.

Gene ID:

TRAT1

Uniprot

Q6PIZ9

Synonyms:

T-cell receptor-associated transmembrane adapter 1 (T-cell receptor-interacting molecule) (TRIM) (pp29/30), TRAT1, TCRIM


Immunogen:

Recombinant Human T-cell receptor-associated transmembrane adapter 1 protein (29-186AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO37162 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human thymus tissue using PACO37162 at dilution of 1:100.