

# BTNL3 Antibody



PACO37174

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

BTNL3

**Uniprot**

Q6UXE8

**Synonyms:**

Butyrophilin-like protein 3 (Butyrophilin-like receptor), BTNL3, BTNLR COLF4100

**Immunogen:**


Recombinant Human Butyrophilin-like protein 3 protein (18-237AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37174 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO37174 at dilution of 1:100.