

PAQR7 Antibody, Biotin conjugated



PACO37245

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Steroid membrane receptor. Binds progesterone in vitro. May be involved in oocyte maturation.

Gene ID:

PAQR7

Uniprot

Q86WK9

Synonyms:

Membrane progesterin receptor alpha (mPR alpha) (Membrane progesterone P4 receptor alpha) (Membrane progesterone receptor alpha) (Progesterone and adipoQ receptor family member 7) (Progesterin and adipoQ receptor family member 7) (Progesterin and adipoQ receptor family member VII), PAQR7, MRPA

Immunogen:

Recombinant Human Membrane progesterin receptor α protein (1-72AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A