

BAGE3 Antibody



PACO37254

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Unknown. Candidate gene encoding tumor antigens.

Gene ID:

BAGE3

Uniprot

Q86Y29

Synonyms:

B melanoma antigen 3 (Cancer/testis antigen 2.3) (CT2.3), BAGE3


Immunogen:

Recombinant Human B melanoma antigen 3 protein (18-109AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung cancer using PACO37254 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO37254 at dilution of 1:100.