

# TRPM8 Antibody, FITC conjugated



PACO37264

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TRPM8

**Uniprot**

Q7Z2W7

**Synonyms:**

Transient receptor potential cation channel subfamily M member 8 (Long transient receptor potential channel 6) (LTrpC-6) (LTrpC6) (Transient receptor potential p8) (Trp-p8), TRPM8, LTRPC6 TRPP8

**Immunogen:**

Recombinant Human Transient receptor potential cation channel subfamily M member 8 protein (1-192AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A