

PACO37294

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Functions as a protein kinase A-anchoring protein that scaffolds the cAMP-dependent protein kinase holoenzyme. May serve as a point of convergence for MAPK and PKA signaling in cilia.

**Gene ID:**

RSPH3

**Uniprot**

Q86UC2

**Synonyms:**

Radial spoke head protein 3 homolog (A-kinase anchor protein RSPH3) (Radial spoke head-like protein 2), RSPH3, RSHL2 RSP3

**Immunogen:**

Recombinant Human Radial spoke head protein 3 homolog protein (1-275AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A