

VASH2 Antibody



PACO37302

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Angiogenesis inhibitor. Inhibits network formation by endothelial cells.

Gene ID:

VASH2

Uniprot

Q86V25

Synonyms:

Vasohibin-2 (Vasohibin-like protein), VASH2, VASHL


Immunogen:

Recombinant Human Vasohibin-2 protein (1-156AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO37302 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO37302 at dilution of 1:100.