

FAM19A1 Antibody



PACO37330

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FAM19A1

Uniprot

Q7Z5A9

Synonyms:

Protein FAM19A1 (Chemokine-like protein Tafa-1), FAM19A1, Tafa1

Immunogen:

Recombinant Human Protein FAM19A1 protein (20-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A