

Product Information

Size:

50ug

Reactivity:

Chimpanzee

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cytokine that binds to LTBR/TNFRSF3. May play a specific role in immune response regulation. Provides the membrane anchor for the attachment of the heterotrimeric complex to the cell surface.

Gene ID:

LTB

Uniprot

Q862Z7

Synonyms:

Lymphotoxin-beta (LT-beta) (Tumor necrosis factor C) (TNF-C) (Tumor necrosis factor ligand superfamily member 3), LTB, TNFC TNFSF3

Immunogen:

Recombinant Chimpanzee Lymphotoxin-β protein (49-244AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A