

PACO37354

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

nucleus, identical protein binding, poly(A) RNA binding.

Gene ID:

RALYL

Uniprot

Q86SE5

Synonyms:

RNA-binding Raly-like protein (hRALYL) (Heterogeneous nuclear ribonucleoprotein C-like 3) (hnRNP core protein C-like 3), RALYL, HNRPCL3


Immunogen:

Recombinant Human RNA-binding Raly-like protein (1-291AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO37354 at dilution of 1:100.